h/p/cosmos airwalk anwenderhandbuch

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Bestellnummer [cos14323]

Autorenteam

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- BODY WEIGHT SUPPORT TREADMILL TRAINING IMPROVES BLOOD GLUCOSE REGULATION IN PERSONS WITH INCOMPLETE SPINAL CORD INJURY. Phillips SM, et al.; J Appl Physiol. 2004 Apr 23 [Epub ahead of print]; PMID: 15107410 [PubMed - as supplied by publisher]
- LOCOMOTOR ACTIVITY IN SPINAL CORD-INJURED PERSONS. Dietz V, et al.; J Appl Physiol. 2004 May;96(5):1954-1960; PMID: 15075315 [PubMed - as supplied by publisher]
- DISTRIBUTED PLASTICITY OF LOCOMOTOR PATTERN GENERATORS IN SPINAL CORD INJURED PATIENTS. Grasso R, et al.; Brain. 2004 May;127(Pt 5):1019-34. Epub 2004 Feb 26; PMID: 14988161 [PubMed - in process]
- TREADMILL TRAINING AND/OR BODY WEIGHT SUPPORT MAY NOT IMPROVE WALKING ABILITY FOLLOWING STROKE. HELBOSTAD. JL.; Aust J Physiother. 2003;49(4):278. No abstract available. PMID: 14714540 [PubMed]
- TREADMILL TRAINING WITH PARTIAL BODY-WEIGHT SUPPORT AFTER TOTAL HIP ARTHROPLASTY: A RANDOMIZED CONTROLLED TRIAL. Hesse S, et al.; Arch Phys Med Rehabil. 2003 Dec;84(12):1767-73. PMID: 14669181 [PubMed - indexed for MEDLINE]
- UPPER AND LOWER EXTREMITY ROBOTIC DEVICES FOR REHABILITATION AND FOR STUDYING MOTOR CONTROL. Hesse S, et al.; Curr Opin Neurol. 2003 Dec;16(6):705-10. PMID: 14624080 [PubMed - in process]
- TREADMILL TRAINING AND BODY WEIGHT SUPPORT FOR WALKING AFTER STROKE. Moseley AM, et al.; Stroke. 2003 Dec;34(12):3006. Epub 2003 Nov 13. Review. No abstract available. PMID: 14615617 [PubMed - indexed for MEDLINE]
- OPTIMAL OUTCOMES OBTAINED WITH BODY-WEIGHT SUPPORT COMBINED WITH TREADMILL TRAINING IN STROKE SUBJECTS. Barbeau H, et al.; Arch Phys Med Rehabil. 2003 Oct;84(10):1458-65. PMID: 14586912 [PubMed - indexed for MEDLINE]
- GAIT RETRAINING POST STROKE. Teasell RW, et al.; Top Stroke Rehabil. 2003 Summer;10(2):34-65. PMID: 13680517 [PubMed - indexed for MEDLINE]
- THE EFFECT OF TREADMILL TRAINING ON MOTOR RECOVERY AFTER A PARTIAL SPINAL CORD COMPRESSION-INJURY IN THE ADULT RAT. Multon S, et al.; J Neurotrauma. 2003 Aug;20(8):699-706. PMID: 12965049 [PubMed - indexed for MEDLINE]
- TREADMILL TRAINING AND BODY WEIGHT SUPPORT FOR WALKING AFTER STROKE. Moseley AM, et al.; Cochrane Database Syst Rev. 2003;(3):CD002840. Review. PMID: 12917932 [PubMed - indexed for MEDLINE]
- TEMPORAL COMPONENTS OF THE MOTOR PATTERNS EXPRESSED BY THE HUMAN SPINAL CORD REFLECT FOOT KINEMATICS. Ivanenko YP, et al.; J Neurophysiol. 2003 Nov;90(5):3555-65. Epub 2003 Jul 09. PMID: 12853436 [PubMed - indexed for MEDLINE]
- THE EFFECT OF TREADMILL TRAINING ON GAIT, BALANCE AND TRUNK CONTROL IN A HEMIPLEGIC SUBJECT: A SINGLE SYSTEM DESIGN. Mudge S, et al.; Disabil Rehabil. 2003 Sep 2;25(17):1000-7. PMID: 12851089 [PubMed - indexed for MEDLINE]
- TEMPOROSPATIAL AND KINEMATIC GAIT ALTERATIONS DURING TREADMILL WALKING WITH BODY WEIGHT SUSPENSION. Threlkeld AJ, et al.; Gait Posture. 2003 Jun;17(3):235-45. PMID: 12770637 [PubMed - indexed for MEDLINE]
- TREADMILL TRAINING WITH PARTIAL BODY WEIGHT SUPPORT AFTER STROKE. Hesse S, et al.; Phys Med Rehabil Clin N Am. 2003 Feb;14(1 Suppl):S111-23. Review. PMID: 12625641 [PubMed - indexed for MEDLINE]
- EXERCISE CAPACITY EARLY AFTER STROKE. Mackay-Lyons MJ, et al.; Arch Phys Med Rehabil. 2002 Dec;83(12):1697-702. PMID: 12474172 [PubMed - indexed for MEDLINE]
- TREADMILL TRAINING WITH PARTIAL BODY WEIGHT SUPPORT AND AN ELECTROMECHANICAL GAIT TRAINER FOR RESTORATION OF GAIT IN SUBACUTE STROKE PATIENTS: A RANDOMIZED CROSSOVER STUDY. Werner C, et al.; Stroke. 2002 Dec;33(12):2895-901. PMID: 12468788 [PubMed - indexed for MEDLINE]
- REHA-STEPPER LOCOMOTION THERAPY IN EARLY REHABILITATION OF PARAPLEGIC PATIENTS
Rupp R, et al.; Biomed Tech (Berl). 2002;47 Suppl 1 Pt 2:708-11. German. PMID: 12465280 [PubMed - indexed for MEDLINE]
- TREADMILL TRAINING WITH PARTIAL BODY WEIGHT SUPPORT AND PHYSIOTHERAPY IN STROKE PATIENTS: A PRELIMINARY COMPARISON. Werner C, et al.; Eur J Neurol. 2002 Nov;9(6):639-44. PMID: 12453080 [PubMed - indexed for MEDLINE]
- RECRUITMENT OF SPINAL MOTOR POOLS DURING VOLUNTARY MOVEMENTS VERSUS STEPPING AFTER HUMAN SPINAL CORD INJURY. Maegele M, et al.; J Neurotrauma. 2002 Oct;19(10):1217-29. PMID: 12427330 [PubMed - indexed for MEDLINE]
- PREMOTOR CORTEX IS INVOLVED IN RESTORATION OF GAIT IN STROKE. Miyai I, et al.; Ann Neurol. 2002 Aug;52(2):188-94. PMID: 12210789 [PubMed - indexed for MEDLINE]

IMPROVED INTRALIMB COORDINATION IN PEOPLE WITH INCOMPLETE SPINAL CORD INJURY FOLLOWING TRAINING WITH BODY WEIGHT SUPPORT AND ELECTRICAL STIMULATION. Field-Fote EC, et al.; Phys Ther. 2002 Jul;82(7):707-15. PMID: 12088467 [PubMed - indexed for MEDLINE]

CONTROL OF FOOT TRAJECTORY IN HUMAN LOCOMOTION: ROLE OF GROUND CONTACT FORCES IN SIMULATED REDUCED GRAVITY. Ivanenko YP, et al.; J Neurophysiol. 2002 Jun;87(6):3070-89. PMID: 12037209 [PubMed - indexed for MEDLINE]

TREADMILL AMBULATION WITH PARTIAL BODY WEIGHT SUPPORT FOR THE TREATMENT OF LOW BACK AND LEG PAIN. Joffe D, et al.; J Orthop Sports Phys Ther. 2002 May;32(5):202-13; discussion 213-5. PMID: 12014824 [PubMed - indexed for MEDLINE]

STEP TRAINING WITH BODY WEIGHT SUPPORT: EFFECT OF TREADMILL SPEED AND PRACTICE PARADIGMS ON POSTSTROKE LOCOMOTOR RECOVERY. Sullivan KJ, et al.; Arch Phys Med Rehabil. 2002 May;83(5):683-91. PMID: 11994808 [PubMed - indexed for MEDLINE]

SUPPORTED TREADMILL TRAINING FOR GAIT AND BALANCE IN A PATIENT WITH PROGRESSIVE SUPRANUCLEAR PALSY. Suteerawattananon M, et al.; Phys Ther. 2002 May; 82(5):485-95. PMID: 11991801 [PubMed - indexed for MEDLINE]

GAIT TRAINING IN HEMIPLEGIA. Mauritz KH.; Eur J Neurol. 2002 May;9 Suppl 1:23-9; discussion 53-61. Review. PMID: 11918646 [PubMed - indexed for MEDLINE]

BODY WEIGHT SUPPORTED TREADMILL AND OVERGROUND TRAINING IN A PATIENT POST CEREBROVASCULAR ACCIDENT. Miller EW. NeuroRehabilitation. 2001;16(3):155-63. PMID: 11790900 [PubMed - indexed for MEDLINE]

LOCOMOTOR THERAPY IN NEUROREHABILITATION. Hesse S.; NeuroRehabilitation. 2001;16(3):133-9. PMID: 11790898 [PubMed - indexed for MEDLINE]

BODY WEIGHT SUPPORT TREADMILL AND OVERGROUND AMBULATION TRAINING FOR TWO PATIENTS WITH CHRONIC DISABILITY SECONDARY TO STROKE. Miller EW, et al.; Phys Ther. 2002 Jan;82(1):53-61. PMID: 11784278 [PubMed - indexed for MEDLINE]

EFFECT OF 15% BODY WEIGHT SUPPORT ON EXERCISE CAPACITY OF ADULTS WITHOUT IMPAIRMENTS. MacKay-Lyons M, et al.; Phys Ther. 2001 Nov;81(11):1790-800. PMID: 11694172 [PubMed - indexed for MEDLINE]

INFLUENCE OF WALKING SPEED ON LOWER LIMB MUSCLE ACTIVITY AND ENERGY CONSUMPTION DURING TREADMILL WALKING OF HEMIPARETIC PATIENTS. Hesse S, et al.; Arch Phys Med Rehabil. 2001 Nov;82(11):1547-50. PMID: 11689974 [PubMed - indexed for MEDLINE]

WALKING TRAINING OF PATIENTS WITH HEMIPARESIS AT AN EARLY STAGE AFTER STROKE: A COMPARISON OF WALKING TRAINING ON A TREADMILL WITH BODY WEIGHT SUPPORT AND WALKING TRAINING ON THE GROUND. Nilsson L, et al.; Clin Rehabil. 2001 Oct;15(5):515-27. PMID: 11594641 [PubMed - indexed for MEDLINE]

ACUTE EFFECTS OF LOCOMOTOR TRAINING ON OVERGROUND WALKING SPEED AND H-REFLEX MODULATION IN INDIVIDUALS WITH INCOMPLETE SPINAL CORD INJURY. Trimble MH, et al.; J Spinal Cord Med. 2001 Summer;24(2):74-80. PMID: 11587422 [PubMed - indexed for MEDLINE]

MODULATION OF LOCOMOTOR-LIKE EMG ACTIVITY IN SUBJECTS WITH COMPLETE AND INCOMPLETE SPINAL CORD INJURY. Dobkin BH, et al.; J Neurol Rehabil. 1995;9(4):183-90. PMID: 11539274 [PubMed - indexed for MEDLINE]

BODY WEIGHT-SUPPORTED TREADMILL TRAINING AFTER STROKE. Hesse S, et al.; Curr Atheroscler Rep. 2001 Jul;3(4):287-94. Review. PMID: 11389793 [PubMed - indexed for MEDLINE]

COMBINED USE OF BODY WEIGHT SUPPORT, FUNCTIONAL ELECTRIC STIMULATION, AND TREADMILL TRAINING TO IMPROVE WALKING ABILITY IN INDIVIDUALS WITH CHRONIC INCOMPLETE SPINAL CORD INJURY. Field-Fote EC.; Arch Phys Med Rehabil. 2001 Jun;82(6):818-24. PMID: 11387589 [PubMed - indexed for MEDLINE]

EQUIPMENT SPECIFICATIONS FOR SUPPORTED TREADMILL AMBULATION TRAINING. Wilson MS, et al.; J Rehabil Res Dev. 2000 Jul-Aug;37(4):415-22. PMID: 11028697 [PubMed - indexed for MEDLINE]

OXYGEN CONSUMPTION DURING TREADMILL WALKING WITH AND WITHOUT BODY WEIGHT SUPPORT IN PATIENTS WITH HEMIPARESIS AFTER STROKE AND IN HEALTHY SUBJECTS. Danielsson A, et al.; Arch Phys Med Rehabil. 2000 Jul;81(7):953-7. PMID: 10896011 [PubMed - indexed for MEDLINE]

TREADMILL TRAINING WITH BODY WEIGHT SUPPORT: ITS EFFECT ON PARKINSON'S DISEASE. Miyai I, et al.; Arch Phys Med Rehabil. 2000 Jul;81(7):849-52. PMID: 10895994 [PubMed - indexed for MEDLINE]

TREADMILL TRAINING WITH PARTIAL BODY WEIGHT SUPPORT IN NONAMBULATORY PATIENTS WITH CEREBRAL PALSY. Schindl MR, et al.; Arch Phys Med Rehabil. 2000 Mar;81(3):301-6. PMID: 10724074 [PubMed - indexed for MEDLINE]

ELECTROMYOGRAPHIC ANALYSIS AND ENERGY EXPENDITURE OF HARNESS SUPPORTED TREADMILL WALKING: IMPLICATIONS FOR KNEE REHABILITATION. Colby SM, et al.; Gait Posture. 1999 Dec;10(3):200-5. PMID: 10567751 [PubMed - indexed for MEDLINE]

LAUFBAND (TREADMILL) THERAPY IN INCOMPLETE PARAPLEGIA AND TETRAPLEGIA. Wernig A, et al.; J Neurotrauma. 1999 Aug;16(8):719-26. PMID: 10511245 [PubMed - indexed for MEDLINE]

- GAIT PATTERN OF SEVERELY DISABLED HEMIPARETIC SUBJECTS ON A NEW CONTROLLED GAIT TRAINER AS COMPARED TO ASSISTED TREADMILL WALKING WITH PARTIAL BODY WEIGHT SUPPORT. Hesse S, et al.; Clin Rehabil. 1999 Oct;13(5):401-10. PMID: 10498347 [PubMed - indexed for MEDLINE]
- THE GAIT OF PATIENTS WITH FULL WEIGHTBEARING CAPACITY AFTER HIP PROSTHESIS IMPLANTATION ON THE TREADMILL WITH PARTIAL BODY WEIGHT SUPPORT, DURING ASSISTED WALKING AND WITHOUT CRUTCHES; Hesse S, et al. Z Orthop Ihre Grenzgeb. 1999 May-Jun;137(3):265-72. German; PMID: 10441834 [PubMed - indexed for MEDLINE]
- TREADMILL WALKING WITH PARTIAL BODY WEIGHT SUPPORT VERSUS FLOOR WALKING IN HEMIPARETIC SUBJECTS. Hesse S, et al.; Arch Phys Med Rehabil. 1999 Apr;80(4):421-7. PMID: 10206604 [PubMed - indexed for MEDLINE]
- A NEW APPROACH TO RETRAIN GAIT IN STROKE PATIENTS THROUGH BODY WEIGHT SUPPORT AND TREADMILL STIMULATION. Visintin M, et al.; Stroke. 1998 Jun;29(6):1122-8. PMID: 9626282 [PubMed - indexed for MEDLINE]
- PARTIAL BODY WEIGHT SUPPORT WITH TREADMILL LOCOMOTION TO IMPROVE GAIT AFTER INCOMPLETE SPINAL CORD INJURY: A SINGLE-SUBJECT EXPERIMENTAL DESIGN. Gardner MB, et al.; Phys Ther. 1998 Apr;78(4):361-74. PMID: 9555919 [PubMed - indexed for MEDLINE]
- HUMAN LUMBOSACRAL SPINAL CORD INTERPRETS LOADING DURING STEPPING. Harkema SJ, et al.; J Neurophysiol. 1997 Feb;77(2):797-811. PMID: 9065851 [PubMed - indexed for MEDLINE]
- MECHANICAL UNWEIGHTING EFFECTS ON TREADMILL EXERCISE AND PAIN IN ELDERLY PEOPLE WITH OSTEOARTHRITIS OF THE KNEE. Mangione KK, et al.; Phys Ther. 1996 Apr;76(4):387-94. PMID: 8606901 [PubMed - indexed for MEDLINE]
- RESTORATION OF GAIT BY COMBINED TREADMILL TRAINING AND MULTICHANNEL ELECTRICAL STIMULATION IN NON-AMBULATORY HEMIPARETIC PATIENTS. Hesse S, et al.; Scand J Rehabil Med. 1995 Dec;27(4):199-204. PMID: 8650503 [PubMed - indexed for MEDLINE]
- A TREADMILL APPARATUS AND HARNESS SUPPORT FOR EVALUATION AND REHABILITATION OF GAIT. Norman KE, et al.; Arch Phys Med Rehabil. 1995 Aug;76(8):772-8. PMID: 7632134 [PubMed - indexed for MEDLINE]
- TREADMILL TRAINING WITH PARTIAL BODY WEIGHT SUPPORT COMPARED WITH PHYSIOTHERAPY IN NONAMBULATORY HEMIPARETIC PATIENTS. Hesse S, et al.; Stroke. 1995 Jun;26(6):976-81. PMID: 7762049 [PubMed - indexed for MEDLINE]
- LAUFBAND THERAPY BASED ON 'RULES OF SPINAL LOCOMOTION' IS EFFECTIVE IN SPINAL CORD INJURED PERSONS. Wernig A, et al.; Eur J Neurosci. 1995 Apr 1;7(4):823-9. Erratum in: Eur J Neurosci 1995 Jun 1;7(6):1429. PMID: 7620630 [PubMed - indexed for MEDLINE]
- RESTORATION OF GAIT IN NONAMBULATORY HEMIPARETIC PATIENTS BY TREADMILL TRAINING WITH PARTIAL BODY-WEIGHT SUPPORT. Hesse S, et al.; Arch Phys Med Rehabil. 1994 Oct;75(10):1087-93. PMID: 7944913 [PubMed - indexed for MEDLINE]
- THE EFFECTS OF PARALLEL BARS, BODY WEIGHT SUPPORT AND SPEED ON THE MODULATION OF THE LOCOMOTOR PATTERN OF SPASTIC PARETIC GAIT. A PRELIMINARY COMMUNICATION. Visintin M, et al.; Paraplegia. 1994 Aug;32(8):540-53. PMID: 7970859 [PubMed - indexed for MEDLINE]
- LAUFBAND LOCOMOTION WITH BODY WEIGHT SUPPORT IMPROVED WALKING IN PERSONS WITH SEVERE SPINAL CORD INJURIES. Wernig A, et al.; Paraplegia. 1992 Apr;30(4):229-38. PMID: 1625890 [PubMed - indexed for MEDLINE]
- INFLUENCE OF BODY WEIGHT SUPPORT ON NORMAL HUMAN GAIT: DEVELOPMENT OF A GAIT RETRAINING STRATEGY. Finch L, et al.; Phys Ther. 1991 Nov;71(11):842-55; discussion 855-6. PMID: 1946621 [PubMed - indexed for MEDLINE]
- THE COMBINED EFFECTS OF CLONIDINE AND CYPROHEPTADINE WITH INTERACTIVE TRAINING ON THE MODULATION OF LOCOMOTION IN SPINAL CORD INJURED SUBJECTS. Fung J, et al.; J Neurol Sci. 1990 Dec;100(1-2):85-93. PMID: 2089144 [PubMed - indexed for MEDLINE]
- THE EFFECTS OF CYPROHEPTADINE ON LOCOMOTION AND ON SPASTICITY IN PATIENTS WITH SPINAL CORD INJURIES. Wainberg M, et al.; J Neurol Neurosurg Psychiatry. 1990 Sep;53(9):754-63. PMID: 2246657 [PubMed - indexed for MEDLINE]
- THE EFFECTS OF BODY WEIGHT SUPPORT ON THE LOCOMOTOR PATTERN OF SPASTIC PARETIC PATIENTS. Visintin M, et al.; Can J Neurol Sci. 1989 Aug;16(3):315-25. PMID: 2766124 [PubMed - indexed for MEDLINE]

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Weight supported treadmill training: Clinical studies [George Chen, Biomedical Engineer, Stanford University, California]

Harness supported treadmill training

This section summarizes the different approaches taken by experimenters applying harness-supported treadmill training to neurologically impaired subjects and the current status of work quantifying the immediate gait response (i.e., gait kinematics, temporal-distance parameters, and symmetry measures) of normal and neurologically impaired subjects to the training technique. Of the three training parameters proposed (i.e., body weight support, treadmill speed, and harness-support compliance), only the effect of treadmill speed on training outcomes has been studied. This study on training speed by Sullivan et al. (2000) is summarized, along with some preliminary work and opinions on harness-support compliance.

Training Approaches

Many approaches have been used to prescribe body weight support, treadmill speed, and manual assistance in the application of harness-supported treadmill training. Hesse et al. (1995) set initial body weight support to 30% based on clinical experience, and reduced it as rapidly as possible to ensure full weight bearing. Treadmill speed was kept deliberately slow to permit longer training sessions and facilitate gait corrections using manual assistance (Hesse et al., 1995). Gardner et al. (1998) set body weight support to the highest level that allowed the subject to achieve heel contact for 10 consecutive steps and maintained this level throughout the study. Treadmill speed was increased in 0.5 mph increments when the subject was able to ambulate without scuffing the paretic foot for 10 consecutive steps during the maximum speed in the prior session (Gardner et al., 1998). Manual assistance was not provided in this study since the subject was able to ambulate independently (Gardner et al., 1998). Visintin et al. (1998) observed subjects walking at 10%, 20%, 30%, and 40% body weight support and selected the percent body weight support that facilitated proper trunk and limb alignment and transfer of weight onto the hemiparetic limb. One or two therapists provided manual assistance, as needed (Visintin et al., 1998). Body weight support was reduced and treadmill speed increased in a stepwise manner as the subject's walking ability improved (Visintin et al., 1998). Seif-Naraghi and Herman (1999) described a number of training approaches, including one that advocated training subjects at the fastest possible speed with as much help from the device and experimenters as needed. Considerable latitude exists in the application of harness-supported treadmill training because subjects are capable of walking at different training parameter settings and with variable amounts of manual assistance during training sessions. Experimenters can choose to reduce body weight support as rapidly as possible or give training speed a higher priority. Manual assistance can be provided only when absolutely necessary (subject cannot practice stepping) or whenever the subject is not producing the "desired" gait kinematics. However, it's often unclear when optimal gait kinematics and, particularly, gait kinetics are produced as training parameters are adjusted.

Gait Kinematics

Gait kinematics have been reported during harness-supported treadmill walking in spastic paretic (Visintin and Barbeau, 1989), hemiparetic (Hesse et al., 1997), and neurologically healthy subjects (Finch et al., 1991). In these studies, joint angular displacements were assessed manually from the monitor screen with the aid of reflective markers (Finch et al., 1991; Visintin and Barbeau, 1989) or qualitatively by raters (Hesse et al., 1997). At 15-60% body weight support, both spastic paretic (Visintin and Barbeau, 1989) and hemiparetic (Hesse et al., 1997) subjects walked more upright with straighter hip and knee alignment during stance. At up to 45% body weight support, the authors noted more normal joint angular displacement profiles (Visintin and Barbeau, 1989) and a more physiologic, plantigrade manner of weight acceptance (Hesse et al., 1997). However, at 60% body weight support, hemiparetic subjects tended to walk on their toes, particularly with the non-paretic limb (Hesse et al., 1997). This amount of body weight support was therefore regarded as unfavorable by the authors (Hesse et al., 1997). Interestingly, body weight support also decreased hip and knee flexion during swing (Hesse et al., 1997). Straighter hip and knee alignment during stance and decreased hip and knee flexion during swing were also observed in neurologically healthy subjects walking with body weight support (Finch et al., 1991). The authors noted that harness constraints limiting the downward excursion of the center of gravity might have contributed to these changes (Finch et al., 1991). Indeed, trunk height was found to increase with body weight support (Finch et al., 1991).

The authors also speculated that a decrease in kinetic energy transfer during terminal stance, resulting in decreased swing momentum and displacement, might have contributed to the decreased hip and knee flexion observed during swing (Finch et al., 1991).

Temporal-Distance and Symmetry Measures

Temporal-distance and symmetry measures during harness-supported treadmill walking have been determined using foot switches (Finch et al., 1991; Hassid et al., 1997; Visintin and Barbeau, 1989) or instrumented overshoe slippers (Hesse et al., 1997; Hesse et al., 1999). In spastic paretic (Visintin and Barbeau, 1989), hemiparetic (Hesse et al., 1997; Hesse et al., 1999), and neurologically healthy subjects (Finch et al., 1991), single limb support time increased (particularly, in the paretic limb(s)) and relative double support time decreased with body weight support. The authors suggested that body weight support might provide a greater stimulus for balance training since single limb support requires the paretic limb to both balance and bear weight while the contralateral limb is in swing (Finch et al., 1991; Hesse et al., 1997; Hesse et al., 1999; Visintin and Barbeau, 1989). Visintin et al. (1989) also reported increased stride length and maximum comfortable walking speed with body weight support. In two studies on hemiparetic subjects by Hesse et al. (1997; 1999), swing symmetry increased with body weight support, but stance, double limb support, and step length symmetry were not affected. However, stance symmetry was greater while walking on the treadmill than walking overground (Hesse et al., 1999). Hassid et al. (1997) reported greater single limb stance symmetry (a measure equivalent to swing symmetry) in hemiparetic subjects while walking on a treadmill as compared to walking overground. Interestingly, this increase in symmetry occurred even at 0% body weight support in some subjects (Hassid et al., 1997). Nevertheless, single limb stance symmetry improved most consistently at 15-30% body weight support and was disrupted at 50% body weight support in some subjects (Hassid et al., 1997). Higher belt speeds did not alter this general improvement in symmetry on the treadmill (Hassid et al., 1997). The authors suggested that treadmill walking might allow hemiparetic subjects to receive more normal and symmetrical step-related sensory feedback, which could enhance locomotion (Hassid et al., 1997).

Harness-Support Compliance

Vertical ground reaction forces and center of mass movement have been measured in spinal cord injured and neurologically healthy subjects as they walked on a treadmill with their weight partially supported by a rigid cable or compliant pneumatic cylinder (Gordon et al., 2000). Support by a rigid cable was found to restrict vertical center of mass movement and reduce peaks in vertical ground reaction force at heel strike and toe off (Gordon et al., 2000). Support by a compliant pneumatic cylinder allowed vertical center of mass movement similar to overground locomotion, but peaks in vertical ground reaction force at heel strike were abnormally greater than at toe off (Gordon et al., 2000). The abnormal vertical ground reaction forces may have been due to viscosity in the pneumatic system since support force still fluctuated with trunk motion. A dynamic regulator was later added to the pneumatic system, which reduced body weight support fluctuation during the gait cycle and allowed vertical ground reaction forces comparable to overground gait adjusted for body weight support (Gordon et al., 2000). The authors concluded that the body weight support system affects center of mass movement and vertical ground reaction forces, which could influence afferent information important to training (Gordon et al., 2000). Wilson et al. (2000) provided equipment specifications on the body weight support system based on clinical feedback and the application of engineering principles. The authors believed a system that does not allow vertical displacement produces an unnatural form of gait, which is not the goal of therapy (Wilson et al., 2000). They recommended a Hooke's law spring system with a stiffness of about 880 N/m, reasoning that it would allow enough vertical motion for normal gait, but not enough to allow the patient to lose posture (Wilson et al., 2000). Moreover, the increase in support when the spring is stretched could provide more feedback and reassurance to patients and a smoother safety catch if they should fall (Wilson et al., 2000).

Training Speed

One study has compared the effectiveness of harness-supported treadmill training using different training speeds. Sullivan et al. (2000) randomly assigned 24 individuals with unilateral stroke to slow (0.5 mph), fast (2.0 mph), and variable (0.5, 1.0, 1.5, 2.0 mph) speed training groups. The subjects received 20 minutes of harness-supported treadmill training for 12 sessions across 4-5 weeks at the prescribed treadmill speed(s) (Sullivan et al., 2000). Even though all groups improved their self-selected overground walking velocity after the training phase, significantly greater improvement occurred in the fast training group (Sullivan et al., 2000). The authors concluded that training speed might be a practice parameter that enhances the effectiveness of harness-supported treadmill training (Sullivan et al., 2000).

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© 2005 h/p/cosmos sports & medical gmbh