

compass 600

smart assist

TECHNICAL DATA

- » For all compass 600 devices



SPECIFICATION

- » Intuitive operator guidance via an 8" touchscreen monitor
- » Quick start/power training without using RFID or chip card
- » RFID/chip card reader
- » Medical LAN Isolator
- » Sinual feedback training

OPTIONS

- » Position sensors for simple, reproducible positioning during testing and training
- » Pulse receiver uncoded

TEST PROCEDURE

- » Range of motion test to determine the physiological amplitude of motion
- » Submaximal testing according to the repetition method, supplemented by an objective assessment of the quality of the test and the subjective assessment of the degree of exertion by the test person using the RPE scale